

# KERN OBF 133

The high-performance compound microscope for every laboratory with fixed, pre-centred Koehler illumination

# KERN®



## Piktogramme

### STANDARD



### OPTION



### Power supply

Input voltage:	100 - 240 V
Mains adapter/adaptor included:	EURO

### Category

Brand:	Optics
--------	--------

### Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	450 x 420 x 275 mm
Gross weight:	7,500 kg
Net weight:	6,500 kg
Shipping method:	Parcel service

### Construction

Dimensions (WxDxH):	394 x 185 x 377 mm
Eye distance (max):	75
Eye distance (min):	50
Optical system:	Finite
Rotation of tube 360°:	yes
Tube:	Trinocular
Tube angle:	30°