

KERN OBS101

The school microscope - For the first steps in microscopy and for use in biology lessons



Construction

Dimensions (WxDxH):	130 x 300 x 310 mm
Optical system:	Finite
Rotation of tube 360°:	yes
Tube:	Monocular
Tube angle:	45°

Piktogramme

STANDARD



Power supply

Input voltage:	100 - 240 V
Mains adapter/adaptor included:	EURO

n/a:	Mains adapter external
------	------------------------

Category

Brand:	Optics
--------	--------

Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	410 x 280 x 210 mm
Gross weight:	3 kg
Net weight:	2,500 kg
Shipping method:	Parcel service
UN number:	3496